





Recombinant Human Oxytocin-neurophysin 1 (OXT), partial

Product Code	CSB-EP017315HUb3
Relevance	Neurophysin 1 specifically binds oxytocin. Oxytocin causes contraction of the smooth muscle of the uterus and of the mammary gland.
Abbreviation	Recombinant Human OXT protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01178
Alias	Ocytocin
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	AAPDLDV RKCLPCGPGGKGRCFGPNICCAEELGCFVGTAEALRCQEENYLPS PCQSGQKACGSGGRCV LGLCCSPDGCHADPACDAEATFSQR
Source	E.coli
Target Names	OXT
Expression Region	32-125aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Mol. Weight	29.6kDa
Protein Length	Partial


Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP017315HUb3 could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) OXT.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP017315HUb3 could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) OXT.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.