



# ACSL1 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA001191LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P33121
<b>Immunogen</b>	Recombinant Human Long-chain-fatty-acid--CoA ligase 1 protein (435-547AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Long-chain-fatty-acid--CoA ligase 1 (EC 6.2.1.3) (Acyl-CoA synthetase 1) (ACS1) (Long-chain acyl-CoA synthetase 1) (LACS 1) (Long-chain acyl-CoA synthetase 2) (LACS 2) (Long-chain fatty acid-CoA ligase 2) (Palmitoyl-CoA ligase 1) (Palmitoyl-CoA ligase 2), ACSL1, FACL1 FACL2 LACS LACS1 LACS2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	ACSL1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.