



ACSL1 Antibody, Biotin conjugated

Product Code	CSB-PA001191LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33121
Immunogen	Recombinant Human Long-chain-fatty-acid--CoA ligase 1 protein (435-547AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Long-chain-fatty-acid--CoA ligase 1 (EC 6.2.1.3) (Acyl-CoA synthetase 1) (ACS1) (Long-chain acyl-CoA synthetase 1) (LACS 1) (Long-chain acyl-CoA synthetase 2) (LACS 2) (Long-chain fatty acid-CoA ligase 2) (Palmitoyl-CoA ligase 1) (Palmitoyl-CoA ligase 2), ACSL1, FACL1 FACL2 LACS LACS1 LACS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACSL1