




CYP1A2 Antibody

Product Code	CSB-PA006396LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05177
Immunogen	Recombinant Human Cytochrome P450 1A2 protein (268-445AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 1A2 (EC 1.14.14.1) (CYPIA2) (Cholesterol 25-hydroxylase) (Cytochrome P(3)450) (Cytochrome P450 4) (Cytochrome P450-P3), CYP1A2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP1A2

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA006396LA01HU at dilution of 1:100