



EPHB3 Antibody, FITC conjugated

Product Code	CSB-PA007731LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54753
Immunogen	Recombinant Human Ephrin type-B receptor 3 protein (261-390AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin type-B receptor 3 (EC 2.7.10.1) (EPH-like tyrosine kinase 2) (EPH-like kinase 2) (Embryonic kinase 2) (EK2) (hEK2) (Tyrosine-protein kinase TYRO6), EPHB3, ETK2 HEK2 TYRO6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EPHB3