




FLI1 Antibody

Product Code	CSB-PA008714LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01543
Immunogen	Recombinant Human Friend leukemia integration 1 transcription factor protein (205-331AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Friend leukemia integration 1 transcription factor (Proto-oncogene Fli-1) (Transcription factor ERGB), FLI1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FLI1

Image



Western Blot

Positive WB detected in: Jurkat whole cell lysate, K562 whole cell lysate


All lanes: FLI1 antibody at 3µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 51, 44, 48, 30 kDa


Observed band size: 51 kDa



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA008714LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA008714LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA008714LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.