



FLI1 Antibody, HRP conjugated

Product Code	CSB-PA008714LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01543
Immunogen	Recombinant Human Friend leukemia integration 1 transcription factor protein (205-331AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Friend leukemia integration 1 transcription factor (Proto-oncogene Fli-1) (Transcription factor ERGB), FLI1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FLI1