



HTR4 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA010895LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q13639 |
| Immunogen | Recombinant Human 5-hydroxytryptamine receptor 4 protein (316-388AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | 5-hydroxytryptamine receptor 4 (5-HT-4) (5-HT4) (Serotonin receptor 4), HTR4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | HTR4 |