



HTR4 Antibody, Biotin conjugated

Product Code	CSB-PA010895LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13639
Immunogen	Recombinant Human 5-hydroxytryptamine receptor 4 protein (316-388AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5-hydroxytryptamine receptor 4 (5-HT-4) (5-HT4) (Serotonin receptor 4), HTR4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HTR4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.