



IPP Antibody, Biotin conjugated

Product Code	CSB-PA011782LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y573
Immunogen	Recombinant Human Actin-binding protein IPP protein (444-515AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-binding protein IPP (Intracisternal A particle-promoted polypeptide) (IPP) (Kelch-like protein 27), IPP, KLHL27
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IPP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.