



IRX4 Antibody, HRP conjugated

Product Code	CSB-PA011836LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78413
Immunogen	Recombinant Human Iroquois-class homeodomain protein IRX-4 protein (287-363AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iroquois-class homeodomain protein IRX-4 (Homeodomain protein IRXA3) (Iroquois homeobox protein 4), IRX4, IRXA3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IRX4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.