




LCT Antibody

Product Code	CSB-PA012827LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09848
Immunogen	Recombinant Human Lactase-phlorizin hydrolase protein (983-1252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lactase-phlorizin hydrolase (Lactase-glycosylceramidase) [Includes: Lactase (EC 3.2.1.108); Phlorizin hydrolase (EC 3.2.1.62)], LCT, LPH
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LCT

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA012827LA01HU at dilution of 1:100