



MAP2K7 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA013416LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O14733 |
| Immunogen | Recombinant Human Dual specificity mitogen-activated protein kinase kinase 7 protein (17-160AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Dual specificity mitogen-activated protein kinase kinase 7 (MAP kinase kinase 7) (MAPKK 7) (EC 2.7.12.2) (JNK-activating kinase 2) (MAPK/ERK kinase 7) (MEK 7) (Stress-activated protein kinase kinase 4) (SAPK kinase 4) (SAPKK-4) (SAPKK4) (c-Jun N-terminal kinase kinase 2) (JNK kinase 2) (JNKK 2), MAP2K7, JNKK2 MEK7 MKK7 PRKMK7 SKK4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | MAP2K7 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |