



MAP3K9 Antibody, Biotin conjugated

Product Code	CSB-PA013434LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80192
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 9 protein (864-1098AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 9 (EC 2.7.11.25) (Mixed lineage kinase 1), MAP3K9, MLK1 PRKE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAP3K9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.