



MAP3K9 Antibody, Biotin conjugated

Product Code	CSB-PA013434LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80192
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 9 protein (864-1098AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 9 (EC 2.7.11.25) (Mixed lineage kinase 1), MAP3K9, MLK1 PRKE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAP3K9