



MC5R Antibody, Biotin conjugated

Product Code	CSB-PA013562LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33032
Immunogen	Recombinant Human Melanocortin receptor 5 protein (1-37AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanocortin receptor 5 (MC5-R) (MC-2), MC5R
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MC5R