



MNX1 Antibody, Biotin conjugated

Product Code	CSB-PA014692LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50219
Immunogen	Recombinant Human Motor neuron and pancreas homeobox protein 1 protein (50-99AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Motor neuron and pancreas homeobox protein 1 (Homeobox protein HB9), MNX1, HLXB9
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MNX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.