





MAP3K11 Antibody

Product Code	CSB-PA01469A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16584
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 11 protein (536-693AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 11 (EC 2.7.11.25) (Mixed lineage kinase 3) (Src-homology 3 domain-containing proline-rich kinase), MAP3K11, MLK3 PTK1 SPRK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K11

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA01469A0Rb at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA01469A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)