



NOTCH1 Antibody, Biotin conjugated

Product Code	CSB-PA015949LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46531
Immunogen	Recombinant Human Neurogenic locus notch homolog protein 1 protein (2428-2555AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenic locus notch homolog protein 1 (Notch 1) (hN1) (Translocation-associated notch protein TAN-1) [Cleaved into: Notch 1 extracellular truncation (NEXT); Notch 1 intracellular domain (NICD)], NOTCH1, TAN1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NOTCH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.