



OLIG1 Antibody, HRP conjugated

Product Code	CSB-PA016328LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAK6
Immunogen	Recombinant Human Oligodendrocyte transcription factor 1 protein (17-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oligodendrocyte transcription factor 1 (Oligo1) (Class B basic helix-loop-helix protein 6) (bHLHb6) (Class E basic helix-loop-helix protein 21) (bHLHe21), OLIG1, BHLHB6 BHLHE21
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OLIG1