



OLIG1 Antibody, FITC conjugated

Product Code	CSB-PA016328LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAK6
Immunogen	Recombinant Human Oligodendrocyte transcription factor 1 protein (17-105AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oligodendrocyte transcription factor 1 (Oligo1) (Class B basic helix-loop-helix protein 6) (bHLHb6) (Class E basic helix-loop-helix protein 21) (bHLHe21), OLIG1, BHLHB6 BHLHE21
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OLIG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.