



PIR Antibody, Biotin conjugated

Product Code	CSB-PA018032LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00625
Immunogen	Recombinant Human Pirin protein (1-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pirin (EC 1.13.11.24) (Probable quercetin 2,3-dioxygenase PIR) (Probable quercetinase), PIR
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PIR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.