



POU3F1 Antibody, FITC conjugated

Product Code	CSB-PA018396LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03052
Immunogen	Recombinant Human POU domain, class 3, transcription factor 1 protein (59-257AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	POU domain, class 3, transcription factor 1 (Octamer-binding protein 6) (Oct-6) (Octamer-binding transcription factor 6) (OTF-6) (POU domain transcription factor SCIP), POU3F1, OCT6 OTF6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	POU3F1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.