



PSMG1 Antibody, Biotin conjugated

Product Code	CSB-PA018921LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95456
Immunogen	Recombinant Human Proteasome assembly chaperone 1 protein (91-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome assembly chaperone 1 (PAC-1) (Chromosome 21 leucine-rich protein) (C21-LRP) (Down syndrome critical region protein 2), PSMG1, C21LRP, DSCR2, PAC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.