



PSMG1 Antibody, Biotin conjugated

Product Code	CSB-PA018921LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95456
Immunogen	Recombinant Human Proteasome assembly chaperone 1 protein (91-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome assembly chaperone 1 (PAC-1) (Chromosome 21 leucine-rich protein) (C21-LRP) (Down syndrome critical region protein 2), PSMG1, C21LRP DSCR2 PAC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMG1