



PTGER4 Antibody, HRP conjugated

Product Code	CSB-PA018973LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35408
Immunogen	Recombinant Human Prostaglandin E2 receptor EP4 subtype protein (354-472AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostaglandin E2 receptor EP4 subtype (PGE receptor EP4 subtype) (PGE2 receptor EP4 subtype) (Prostanoid EP4 receptor), PTGER4, PTGER2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTGER4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.