



PTPRF Antibody, HRP conjugated

Product Code	CSB-PA019053LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10586
Immunogen	Recombinant Human Receptor-type tyrosine-protein phosphatase F protein (688-970AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-type tyrosine-protein phosphatase F (EC 3.1.3.48) (Leukocyte common antigen related) (LAR), PTPRF, LAR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PTPRF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.