



PTPRF Antibody, FITC conjugated

Product Code	CSB-PA019053LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10586
Immunogen	Recombinant Human Receptor-type tyrosine-protein phosphatase F protein (688-970AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-type tyrosine-protein phosphatase F (EC 3.1.3.48) (Leukocyte common antigen related) (LAR), PTPRF, LAR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PTPRF