



STX1A Antibody, HRP conjugated

Product Code	CSB-PA022895LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16623
Immunogen	Recombinant Human Syntaxin-1A protein (1-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-1A (Neuron-specific antigen HPC-1), STX1A, STX1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	STX1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.