



SYK Antibody, Biotin conjugated

Product Code	CSB-PA022998LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43405
Immunogen	Recombinant Human Tyrosine-protein kinase SYK protein (260-361AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase SYK (EC 2.7.10.2) (Spleen tyrosine kinase) (p72-Syk), SYK
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SYK