



TFE3 Antibody, FITC conjugated

Product Code	CSB-PA023428LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19532
Immunogen	Recombinant Human Transcription factor E3 protein (1-110AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor E3 (Class E basic helix-loop-helix protein 33) (bHLHe33), TFE3, BHLHE33
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TFE3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.