



TPT1 Antibody, Biotin conjugated

Product Code	CSB-PA024134LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13693
Immunogen	Recombinant Human Translationally-controlled tumor protein (1-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Translationally-controlled tumor protein (TCTP) (Fortilin) (Histamine-releasing factor) (HRF) (p23), TPT1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TPT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.