



SLC6A8 Antibody, Biotin conjugated

Product Code	CSB-PA04679D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48029
Immunogen	Recombinant Human Sodium- and chloride-dependent creatine transporter 1 protein (582-635AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium- and chloride-dependent creatine transporter 1 (CT1) (Creatine transporter 1) (Solute carrier family 6 member 8), SLC6A8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC6A8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.