




API5 Antibody


Product Code	CSB-PA05559A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZZ5
Immunogen	Recombinant Human Apoptosis inhibitor 5 protein (396-524AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptosis inhibitor 5 (API-5) (Antiapoptosis clone 11 protein) (AAC-11) (Cell migration-inducing gene 8 protein) (Fibroblast growth factor 2-interacting factor) (FIF) (Protein XAGL), API5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	API5

Image



Western Blot

Positive WB detected in: HL60 whole cell lysate, K562 whole cell lysate, Mouse spleen tissue
 All lanes: API5 antibody at 3µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 60, 50, 57, 51, 58 kDa
 Observed band size: 60, 50 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA05559A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.