



API5 Antibody, Biotin conjugated

Product Code	CSB-PA05559D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZZ5
Immunogen	Recombinant Human Apoptosis inhibitor 5 protein (396-524AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptosis inhibitor 5 (API-5) (Antiapoptosis clone 11 protein) (AAC-11) (Cell migration-inducing gene 8 protein) (Fibroblast growth factor 2-interacting factor) (FIF) (Protein XAGL), API5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	API5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.