



HINT1 Antibody, HRP conjugated

Product Code	CSB-PA05639B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49773
Immunogen	Recombinant Human Histidine triad nucleotide-binding protein 1 protein (2-126AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histidine triad nucleotide-binding protein 1 (EC 3.-.-) (Adenosine 5'-monophosphoramidase) (Protein kinase C inhibitor 1) (Protein kinase C-interacting protein 1) (PKCI-1), HINT1, HINT PKCI1 PRKCNH1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HINT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.