



NFKB2 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA10449B0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q00653 |
| Immunogen | Recombinant Human Nuclear factor NF-kappa-B p100 subunit protein (263-439AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Nuclear factor NF-kappa-B p100 subunit (DNA-binding factor KBF2) (H2TF1) (Lymphocyte translocation chromosome 10 protein) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2) (Oncogene Lym-10) (Lyt10) [Cleaved into: Nuclear factor NF-kappa-B p52 subunit], NFKB2, LYT10 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | NFKB2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |