



NFKB2 Antibody, FITC conjugated

Product Code	CSB-PA10449C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00653
Immunogen	Recombinant Human Nuclear factor NF-kappa-B p100 subunit protein (263-439AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear factor NF-kappa-B p100 subunit (DNA-binding factor KBF2) (H2TF1) (Lymphocyte translocation chromosome 10 protein) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2) (Oncogene Lym-10) (Lyt10) [Cleaved into: Nuclear factor NF-kappa-B p52 subunit], NFKB2, LYT10
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NFKB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.