



APOH Antibody, Biotin conjugated

Product Code	CSB-PA11079D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02749
Immunogen	Recombinant Human Beta-2-glycoprotein 1 protein (119-345AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-2-glycoprotein 1 (APC inhibitor) (Activated protein C-binding protein) (Anticardiolipin cofactor) (Apolipoprotein H) (Apo-H) (Beta-2-glycoprotein I) (B2GPI) (Beta(2)GPI), APOH, B2G1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	APOH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.