




MMP11 Antibody

Product Code	CSB-PA12369A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24347
Immunogen	Recombinant Human Stromelysin-3 protein (198-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stromelysin-3 (SL-3) (ST3) (EC 3.4.24.-) (Matrix metalloproteinase-11) (MMP-11), MMP11, STMY3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP11

Image



Immunofluorescent analysis of A549 cells using CSB-PA12369A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)