



IL12A Antibody, Biotin conjugated

Product Code	CSB-PA16819D0Rb
Abbreviation	Interleukin-12 subunit alpha
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29459
Immunogen	Recombinant Human Interleukin-12 subunit alpha protein (23-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Cytokine that can act as a growth factor for activated T and NK cells, enhance the lytic activity of NK/lymphokine-activated Killer cells, and stimulate the production of IFN-gamma by resting PBMC.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-12 subunit alpha (IL-12A) (Cytotoxic lymphocyte maturation factor 35 kDa subunit) (CLMF p35) (IL-12 subunit p35) (NK cell stimulatory factor chain 1) (NKSF1), IL12A, NKSF1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL12A