



IL12A Antibody, Biotin conjugated

Product Code	CSB-PA16819D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29459
Immunogen	Recombinant Human Interleukin-12 subunit alpha protein (23-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-12 subunit alpha (IL-12A) (Cytotoxic lymphocyte maturation factor 35 kDa subunit) (CLMF p35) (IL-12 subunit p35) (NK cell stimulatory factor chain 1) (NKSF1), IL12A, NKSF1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL12A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.