



PDE4B Antibody, Biotin conjugated

Product Code	CSB-PA21029D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07343
Immunogen	Recombinant Human cAMP-specific 3',5'-cyclic phosphodiesterase 4B protein (27-215AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cAMP-specific 3',5'-cyclic phosphodiesterase 4B (EC 3.1.4.53) (DPDE4) (PDE32), PDE4B, DPDE4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDE4B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.