



ROR2 Antibody, HRP conjugated

Product Code	CSB-PA21119B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01974
Immunogen	Recombinant Human Tyrosine-protein kinase transmembrane receptor ROR2 protein (752-926AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase transmembrane receptor ROR2 (EC 2.7.10.1) (Neurotrophic tyrosine kinase, receptor-related 2), ROR2, NTRKR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ROR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.