




CNP Antibody

Product Code	CSB-PA21939A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09543
Immunogen	Recombinant Human 2',3'-cyclic-nucleotide 3'-phosphodiesterase protein (226-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2',3'-cyclic-nucleotide 3'-phosphodiesterase (CNP) (CNPase) (EC 3.1.4.37), CNP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CNP

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA21939A0Rb at dilution of 1:100