



CNP Antibody, FITC conjugated

Product Code	CSB-PA21939C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09543
Immunogen	Recombinant Human 2',3'-cyclic-nucleotide 3'-phosphodiesterase protein (226-415AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2',3'-cyclic-nucleotide 3'-phosphodiesterase (CNP) (CNPase) (EC 3.1.4.37), CNP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CNP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.