




POLR2M Antibody

Product Code	CSB-PA316685LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6EEV4
Immunogen	Recombinant Human DNA-directed RNA polymerase II subunit GRINL1A, isoforms 4/5 protein (2-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-directed RNA polymerase II subunit GRINL1A, isoforms 4/5 (DNA-directed RNA polymerase II subunit M, isoforms 4/5), POLR2M, GRINL1A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	POLR2M

Image



Immunofluorescent analysis of A549 cells using CSB-PA316685LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)