



Invasin Antibody, Biotin conjugated

Product Code	CSB-PA321672LD01YAQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19196
Immunogen	Recombinant Yersinia enterocolitica Invasin protein (651-835AA)
Raised In	Rabbit
Species Reactivity	Yersinia enterocolitica
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Invasin
Immunogen Species	Yersinia enterocolitica
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.