



pbpA Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA401992HA01CWV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A5I6G4 |
| Immunogen | Recombinant Clostridium botulinum Penicillin-binding protein 1A protein (663-830AA) |
| Raised In | Rabbit |
| Species Reactivity | Clostridium botulinum |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Penicillin-binding protein 1A (PBP1a) [Includes: Penicillin-insensitive transglycosylase (EC 2.4.1.129) (Peptidoglycan TGase); Penicillin-sensitive transpeptidase (EC 3.4.16.4) (DD-transpeptidase)], pbpA |
| Immunogen Species | Clostridium botulinum |
| Research Area | Others |
| Target Names | pbpA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |