




TBX2 Antibody

Product Code	CSB-PA613386LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13207
Immunogen	Recombinant Human T-box transcription factor TBX2 protein (299-573AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-box transcription factor TBX2 (T-box protein 2), TBX2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TBX2

Image



Western Blot
 Positive WB detected in: Rat lung tissue
 All lanes: TBX2 antibody at 2.5µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 76 kDa
 Observed band size: 76 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.